

Ever *Sheer*

SIGVARIS
LIFE FOR LEGS



20-30 mmHg available
in these styles and colors:

- Calf-length
- Thigh-high with Grip-top
- Pantyhose
- Open toe (Calf & Thigh-high) – **New Style!**
Open toe available in Natural and Suntan only

- Dark Navy – **New Color!**
- Natural
- Suntan
- Mocha – **New Color!**
- Black

**Where efficacy
meets fashion!**

**Where efficacy
meets fashion**

SIGVARIS, Inc.
1119 Highway 74 South
Peachtree City, GA
30269 USA
1.800.322.7744

SIGVARIS Corp.
4535 Dobrin
Ville St-Laurent, Quebec
H4R 2L8 Canada
1.800.363.4999

For a certified SIGVARIS dealer near you, please visit
www.sigvarisusa.com www.sigvaris.ca

SIGVARIS, LIFE FOR LEGS, and PFS are registered trademarks of Ganzoni & Cie AG, CH 9014 St.Gallen/
Switzerland, in many countries worldwide. © Copyright 2009 SIGVARIS, Inc. *PCT/FR2009/0560760

LIT-780BROD 12.09

Ever *Sheer* Medical Compression Hosiery

SIGVARIS EverSheer Medical Compression Hosiery for Women features double-covered yarns and a special hi-tech knit construction to provide a best-in-class sheer and fashionable look. EverSheer combines all-day softness and comfort with greatest possible durability!

EverSheer Features & Benefits

- **The only sheer 20-30 mmHg compression hosiery with double-covered spandex (lay-in):**
 - **No bare spandex to irritate the skin**
 - **Easier to don and remove**
- **Excellent breathability and moisture management**
- **Feels soft on the skin**
- **New design using our unique PFS™ System to ensure a comfortable fit with no binding behind the knee**

Superior Thigh Band

Benefits of patented* Sensinnov thigh band with silicone Grip-top:

- **Stays in place with less skin irritation**
- **Excellent moisture management**
- **Comfortable for everyday wear**



Fashionable, yet clinically efficient

Graduated compression stockings are universally recommended for both surgical and non-surgical treatment of varicose, telangiectasias and reticular leg veins. As you know, patients often experience hyperpigmentation, bruising, swelling, and discomfort following vein procedures. Furthermore, many patients are concerned about fashion. Let SIGVARIS help you in satisfying your fashion-conscious patient with our new 20-30 mmHg EverSheer medical compression hosiery for Women compression stockings.

Keep in mind, 76% of patients wearing compression had a satisfying outcome in vessel disappearance as compared to 56.5% of the patients without compression.²

Did you know...

Recent clinical studies have proven that wearing 20-30 mmHg graduated compression stockings for three weeks post sclerotherapy significantly reduces hyperpigmentation and improves overall vessel appearance.^{1 2 3}

EverSheer

782 EverSheer Key Features

| |
|--|
| Sheerness |
| <i>The only compression hosiery with consistent sheerness from ankle through thigh as proven by monochromic picture tests</i> |
| Softness |
| <i>EverSheer is the only 20-30 mmHg compression hosiery with double-covered lay-in yarns, providing a softer feel with less irritation</i> |
| Top Band (Sensinnov) |
| <i>Patented breathable thigh band stays in place with less skin irritation during wear and removal</i> |
| Moisture Management |
| <i>EverSheer enjoys excellent breathability and moisture management characteristics, allowing moisture to wick away from the skin</i> |
| Ease of Donning |
| <i>Double-covered lay-in yarns make EverSheer easier to don because bare spandex never touches the skin</i> |
| Sizing Assortment |
| <i>EverSheer and the Precise Fit Sizing™ System provide ready-to-wear sizes that fit 97% of female population</i> |

Where Efficacy Meets Fashion

1. Efficacy of Graduated Compression Stockings for an Additional 3 Weeks after Sclerotherapy Treatment of Reticular and Telangiectatic Leg Veins. Noothei, Cadag, Magpanty, and Goldman. Dermatol Surg 2009; 35:53-58
 2. Compression after sclerotherapy for telangiectasias and reticular leg veins: A randomized controlled study. Kern, Ramelet, Wütschert, and Hayoz. Journal of Vascular Surgery, 2007; 45:1212-16
 3. Post-Sclerotherapy Compression: Controlled comparative Study of Duration of Compression and Its Effects on Clinical Outcome. Weiss, Sadick, Goldman, Weiss. Dermatol Surg 1999; 25:105-108